



## Stencil, Squeegee & Misprint Cleaner ScreenClean™ Presaturated Wipes





- Engineered for removing inks
- ESD-safe, residue-free
- Plastic-safe
- Excellent replacement for methylene, chloride, acetone & other ink-removing solvents
- No rip. No lint. Heavy duty wipe.
- Unique “slam-shut” package for extra long service life
- Unlimited shelf-life

### Product Information

Stencil, Squeegee & Misprint Cleaner - *ScreenClean* is a cost-effective cleaner for removing thick-film pastes, partially cured epoxies and most inks used in screen printing operations. Using a strong, shred-proof disposable wipe saturated in a powerful cleaning fluid, it quickly breaks down the molecular bonds in inks and pastes for fast and enhanced cleaning. It is safe on all types of screens and plastics and will not cause damaging fragile emulsions. Its slow drying time allows for a more thorough cleaning. Use for perfectly clean stencil screens and associated machinery, squeegees and work surfaces.

*AVAILABILITY: Stocked by MicroCare distributors world-wide. Contact MicroCare for a distributor near you. For more information contact MicroCare directly or visit [www.MicroCare.com](http://www.MicroCare.com).*

Technical Details and Cleaning Data	
Boiling Point	120 °C / 249 °F
Chemical Family	Hydrocarbon Solvent
Cleaning Strength (KB value)	115 (est.)
Evaporation Rate	Slow
Odor	Mild Alcohol
Percent Volatile	100%
Health, Safety and Environmental Data	
Flammable	Yes
Flashpoint	48.9 °C / 120 °F
NFPA	Health: 1 Fire: 2 Instability: 0
RoHS and WEEE Compliant	Yes
REACH	Yes*
US EPA SNAP-Approved	Yes
GHS / HCS Pictogram	
Packaging	
Minicube <sup>2</sup> - 6 lb / 2.72 kg	MCC-CDIG
Cube <sup>2</sup> - 32 lb / 14.5 kg	MCC-CDIP
Drum <sup>1</sup> - 438 lb / 198.6 kg	MCC-CDID
Presaturated Wipes <sup>3</sup> 443 mls. fluid content	100 - 5" x 8" Wipes MCC-CDIW
Recycle	

\*Contact MicroCare for specific details

## Discover Perfectly Clean [www.MicroCare.com](http://www.MicroCare.com)

### MicroCare Corporation

595 John Downey Drive  
New Britain, CT 06051 USA  
CAGE: OATV9  
Tel: +1 860 827 0626  
Toll Free: 1 800 638 0125  
Email: TechSupport@MicroCare.com

### MicroCare Europe Bvba

Vekestraat 29 B11  
Industriezone 't Sas  
1910 Kampenhout, Belgium  
Tel: +32 2 251 95 05  
Fax: +32 2 400 96 39  
Email: EuroSales@MicroCare.com

### MicroCare Asia Pte Ltd

102E Pasir Panjang Road  
#03-01 Citilink Warehouse Complex  
Singapore 118529  
Tel: +65 6271 0182  
Email: TechSupport@MicroCare.sg